

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-FeLV (Clone #8D11)
CATALOG NUMBER:	MFLV-20M-8D11
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% Sodium Azide
STORAGE:	Store between 2 and 8 °C for up to 30 days. For long-term storage, aliquot and store between -20 °C and -80 °C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with feline leukemia virus (FeLV) p27 as determined by ELISA, LFD and western blot. Suggest using this as ELISA detection antibody when capturing with MFLV-20M-CM1. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR CAT OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

MFLV-20M-8D11_Rev060121