



## **PRODUCT PROFILE**

- PRODUCT:** Affinity Purified Chicken anti-DYKDDDDK
- CATALOG NUMBER:** CFLG-45A
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% sodium azide
- STORAGE:** 2 - 8° C
- EXPIRATION:** 24 months after shipping date.
- PROCEDURE:** Chickens were immunized with the DYKDDDDK peptide conjugated to KLH. Eggs were collected and the IgY fraction prepared. Antibody was immunoaffinity purified using peptide coupled to a different carrier. Antibody concentration was determined using an absorbance at 290 nm: 1.32 equals 1.0 mg of IgY.
- SPECIFICATIONS:** Optimal working dilution should be determined for your particular assay conditions.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***