

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Goat anti-Human Fibronectin - Antiserum

CATALOG NUMBER: GFCT-80

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Goats were immunized with highly purified human fibronectin and the resulting antiserum was collected. The antiserum was delipidated, solid phase adsorbed and 0.45 um filtered. Total protein measurement was determined by refractometer. Specific antibody content is determined immunoaffinity purification off an antigen containing immunosorbent.

SPECIFICATIONS: This antiserum reacts with human fibronectin as determined by IEP and Ouchterlony techniques. This antiserum may cross react with fibronectin from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

GFCT-80_Rev030119